



Atorolimumab Biosimilar - Anti-RHD, CD240D mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1220
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG3-kappa
SOURCE	CAS 202833-08-7

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Atorolimumab,P3x22914G4,RHD, CD240D,anti-RHD, CD240D

