



Glebatumumab Biosimilar - Anti-GPNMB extracellular domain mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1238
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2-kappa
SOURCE	CAS 1020264-78-1

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Glebatumumab,CDX-011 (DOX),CR011,glebatumumab vedotin (ADC),GPNMB extracellular domain,anti-GPNMB extracellular domain

