



Indatuximab Biosimilar - Anti-SDC1, CD138 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1277
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	CAS 1238517-16-2

PRODUCT INFORMATION

ORIGIN SPECIES	Chimeric
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Indatuximab,BT-062,nBT062-DM4,SDC1, CD138,anti-SDC1, CD138

