



Caplacizumab Biosimilar - Anti-VWF A1 domain mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1278
CLONALITY	Monoclonal Antibody
ISOTYPE	VH-VH
SOURCE	CAS 915810-67-2

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Caplacizumab,ALX-0081,PMP12A2h1-linker AAA-PMP12A2h1,caplacizumab-yhdp,VWF A1 domain,anti-VWF A1 domain

