



Actoxumab Biosimilar - Anti-toxin A mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1285
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS 1245634-25-6

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Actoxumab,CDA 1,MBL-CDA1,MDX-066,MK-3415A (combination of actoxumab and bezlotoxumab),toxin A,anti-toxin A

