



Bezlotoxumab Biosimilar - Anti-toxin B mAb - Research Grade

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PX-TA1286 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1-kappa |
| SOURCE | CAS 1246264-45-8 |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | Homo sapiens |
| EXPRESSION SYSTEM | Mammalian cells |
| BUFFER | PBS buffer PH7.5 |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | Bezlotoxumab,MK-6072,CDB-1,MDX-1388,MK-3415A (combination of actoxumab and bezlotoxumab),formerly CDB 1,toxin B,anti-toxin B |

