



Placulumab Biosimilar - Anti-TNFSF2, TNF-alpha, TNFA mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1302
CLONALITY	Monoclonal Antibody
ISOTYPE	[V-kappa]2-Fc
SOURCE	CAS 945781-29-3

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Placulumab,ART621,CEP-37247,PN0621,TNFSF2, TNF-alpha, TNFA,anti-TNFSF2, TNF-alpha, TNFA

