



Vorsetuzumab Biosimilar - Anti-CD70, TNFSF7 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1304
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS 1165740-62-4

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Vorsetuzumab,SGN-70,h1F6,CD70, TNFSF7,anti-CD70, TNFSF7

