



Pritoxaximab Biosimilar - Anti-shiga toxin type 1 B subunit mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1318
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS 1351470-16-0

PRODUCT INFORMATION

ORIGIN SPECIES	Chimeric
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Pritoxaximab,Shigamabs-(drug name for the combination of pritoxaximab and setoxaximab),caStx1,shiga toxin type 1 B subunit,anti-shiga toxin type 1 B subunit

