



Enfortumab Biosimilar - Anti-PVRL4 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1333
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS 1346452-25-2

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Enfortumab,AGS-22CE,AGS-22M,AGS-22M6E,unconjugated:AGS-22C3 or AGSM6,PVRL4,anti-PVRL4

