



Lifastuzumab Biosimilar - Anti-SLC34A2 mAb - Research Grade

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PX-TA1334 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1-kappa |
| SOURCE | CAS 1401812-88-1 |

PRODUCT INFORMATION

| | |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ORIGIN SPECIES | Humanized |
| EXPRESSION SYSTEM | Mammalian cells |
| BUFFER | PBS buffer PH7.5 |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | Lifastuzumab, DNIB0600A (conjugate), MNIB2126A (non conjugate), SLC34A2, anti-SLC34A2 |

