



Pasotuxizumab Biosimilar - Anti-CD3E, FOLH1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1353
CLONALITY	Monoclonal Antibody
ISOTYPE	(scFv-heavy-kappa)-(scFv-heavy-lambda) His-Tag
SOURCE	CAS 1442657-12-6

PRODUCT INFORMATION

ORIGIN SPECIES	Chimeric,Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Pasotuxizumab,BAY 2010112,BITE MT112,MT-112,CD3E, FOLH1,anti-CD3E, FOLH1

