



## Vanucizumab Biosimilar - Anti-ANGPT2, VEGFA, mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1355
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	VH-C-kappa-CH2-CH3_V-lambda-CH1_H-GAMMA-1_L-kappa
<b>SOURCE</b>	CAS 1448221-05-3

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Humanized
<b>EXPRESSION SYSTEM</b>	Mammalian cells
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Vanucizumab,RG-7221,RG7221,RO5520985,ANGPT2, VEGFA, ,anti-ANGPT2, VEGFA,

