



Brolucizumab Biosimilar - Anti-VEGFA mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1374
CLONALITY	Monoclonal Antibody
ISOTYPE	scFv-kappa-heavy
SOURCE	CAS 1531589-13-5

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Brolucizumab,ESBA-1008,ESBA1008,RTH258,VEGFA,anti-VEGFA

