



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Monalizumab Biosimilar - Anti-KLRC1, NKG2A, CD159a, CD94 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1392
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	CAS 1228763-95-8

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Monalizumab,IPH-2201,NN-8765,NNC-0141-0000-0100,NNC-141-01000,anti-NKG2A,humZ270,KLRC1, NKG2A, CD159a, CD94,anti-KLRC1, NKG2A, CD159a, CD94



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.