



Emicizumab Biosimilar - Anti-F10, F9 activated form mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1394
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	CAS 1610943-06-0

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Emicizumab,ACE-910,RG6013,hBS910,F10, F9 activated form,anti-F10, F9 activated form

