



Prezalumab Biosimilar - Anti-ICOSLG, B7-H2, B7RP-1, CD275 mAb - Research Grade

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PX-TA1417 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG2-kappa |
| SOURCE | CAS 1523164-68-2 |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | Homo sapiens |
| EXPRESSION SYSTEM | Mammalian cells |
| BUFFER | PBS buffer PH7.5 |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | Prezalumab,AMG 557,MEDI-5872,ICOSLG, B7-H2, B7RP-1, CD275,anti-ICOSLG, B7-H2, B7RP-1, CD275 |

