



Suptavumab Biosimilar - Anti-RSV glycoprotein F mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1429
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS 1629615-23-1

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Suptavumab,H1H3592P3,REGN 2222,SAR-438584,SAR438584,RSV glycoprotein F,anti-RSV glycoprotein F

