



Duvortuxizumab Biosimilar - Anti-CD19, CD3E mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1432
CLONALITY	Monoclonal Antibody
ISOTYPE	scFv -h-CH2-CH3(_scFv)_h-CH2-CH3
SOURCE	CAS 1831098-91-9

PRODUCT INFORMATION

ORIGIN SPECIES	Chimeric,Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Duvortuxizumab,JNJ-64052781,MGD011,RES192M1.2,CD19, CD3E,anti-CD19, CD3E

