



## Lutikizumab Biosimilar - Anti-IL1A , IL1B mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1440
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	[VH-VH'-H-Gamma1-VL-VL'-C-kappa]-dimer
<b>SOURCE</b>	CAS 1791411-57-8

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Humanized
<b>EXPRESSION SYSTEM</b>	Mammalian cells
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Lutikizumab,ABT-981,IL1A , IL1B,anti-IL1A , IL1B

