



## Flotetuzumab Biosimilar - Anti-CD3E, IL3RA, CD123 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1456
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	VH-VL-VH'-VL'
<b>SOURCE</b>	CAS 1664355-28-5

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Humanized
<b>EXPRESSION SYSTEM</b>	Mammalian cells
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Flotetuzumab, MGD-006, MGD006, RES234, S80880, CD3E, IL3RA, CD123, anti-CD3E, IL3RA, CD123

