



Selicrelumab Biosimilar - Anti-CD40, TNFRSF5 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1465
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2-kappa
SOURCE	CAS 1622140-49-1

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Selicrelumab,CP-870,893,RG-7876,RO-7009789,CD40, TNFRSF5,anti-CD40, TNFRSF5

