



Gimsilumab Biosimilar - Anti-CSF2 , GM-CSF mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1478
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS 1648796-29-5

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Gimsilumab,MORAb-022,CSF2 , GM-CSF,anti-CSF2 , GM-CSF

