



Mosunetuzumab Biosimilar - Anti-CD3E, MS4A1, CD20 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1482
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS 1905409-39-3

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Mosunetuzumab,BTCT4465A,RG-7828,RO7030816,CD3E, MS4A1, CD20,anti-CD3E, MS4A1, CD20

