

## Etigilimab Biosimilar - Anti-TIGIT, VSTM3 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1507
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1-kappa
<b>SOURCE</b>	CAS 2044984-83-8

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Humanized
<b>EXPRESSION SYSTEM</b>	Mammalian cells
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	<p>Store at 4 °C for short term. Freeze and store at -80 °C for long term.</p> <p>Addition of 20 to 40 % glycerol may improve cryoprotection.</p> <p>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.</p> <p>Aliquot the product according to your needs to avoid freeze/thaw cycles.</p>
<b>ALIAS/SYNONYMS</b>	Etigilimab,OMP-313M32,TIGIT, VSTM3,anti-TIGIT, VSTM3

