



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Cemiplimab Biosimilar - Anti-PDCD1, PD1, CD279 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1524
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG4-kappa
<b>SOURCE</b>	CAS 1801342-60-8

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Homo sapiens
<b>EXPRESSION SYSTEM</b>	Mammalian cells
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Cemiplimab, REGN-2810, SAR-439684, cemiplimab-rwlc, PDCD1, PD1, CD279, anti-PDCD1, PD1, CD279



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**