



Gancotamab Biosimilar - Anti-ERBB2, EGFR2, CD340 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1542
CLONALITY	Monoclonal Antibody
ISOTYPE	scFv-heavy-lambda
SOURCE	CAS 1509928-00-0

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Gancotamab,MM-302,ERBB2, EGFR2, CD340,anti-ERBB2, EGFR2, CD340

