



Teclistamab Biosimilar - Anti-BCMA-CD3 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1550
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-lambda
SOURCE	CAS 2119595-80-9

PRODUCT INFORMATION

ORIGIN SPECIES	Chimeric, Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Teclistamab ,JNJ-64007957,MS4A1,anti-MS4A1

