



Plamotamab Biosimilar - Anti-CD3E, MS4A1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1563
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa / scFv-h-CH2-CH3
SOURCE	CAS 2138442-31-4

PRODUCT INFORMATION

ORIGIN SPECIES	Chimeric,Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Plamotamab ,XmAb-13676,CD3E, MS4A1,anti-CD3E, MS4A1

