



Tidutamab Biosimilar - Anti-CD3 epsilon, SSTR2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1565
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa / scFv-h-CH2-CH3
SOURCE	CAS 2148354-90-7

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Tidutamab ,XmAb-18087,CD3 epsilon, SSTR2,anti-CD3 epsilon, SSTR2

