



Avdoralimab Biosimilar - Anti-C5AR1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1578
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS 2226393-85-5

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Avdoralimab ,IPH-5401,NN-8210,NNC-0215-0384,C5AR1,anti-C5AR1

