



Gremubamab Biosimilar - Anti-PcrV mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1591
CLONALITY	Monoclonal Antibody
ISOTYPE	VH-CH1-VH-V-KAPPA-CH2-CH3
SOURCE	CAS 1800381-36-5

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens, Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Gremubamab ,MEDI3902,PcrV,anti-PcrV

