



Manelimab Biosimilar - Anti-CD274, PD-L1, B7-H1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1594
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-lambda
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Manelimab ,BCD-135,CD274, PD-L1, B7-H1,anti-CD274, PD-L1, B7-H1

