



## Sonelokimab Biosimilar - Anti-ALB,IL17A mAb - Research Grade

### OVERVIEW

---

|                          |                     |
|--------------------------|---------------------|
| <b>PRODUCT REFERENCE</b> | PX-TA1606           |
| <b>CLONALITY</b>         | Monoclonal Antibody |
| <b>ISOTYPE</b>           | VH-VH'-VH           |
| <b>SOURCE</b>            |                     |

### PRODUCT INFORMATION

---

|                           |   |
|---------------------------|---|
| <b>ORIGIN SPECIES</b>     | Humanized   |
| <b>EXPRESSION SYSTEM</b>  | Mammalian cells   |
| <b>BUFFER</b>             | PBS buffer PH7.5  |
| <b>STORAGE CONDITIONS</b> | Store at 4 °C for short term. Freeze and store at -80 °C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| <b>ALIAS/SYNONYMS</b>     | Sonelokimab ,M-1095,MSB-0010841,ALB,IL17A,anti-ALB,IL17A  |

