



## Teplizumab Biosimilar - Anti-CD3 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1617
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	
<b>SOURCE</b>	CAS 876387-05-2

### PRODUCT INFORMATION

---

#### ORIGIN SPECIES

**EXPRESSION SYSTEM** Mammalian cells

**BUFFER** PBS buffer PH7.5

**STORAGE CONDITIONS** Store at 4 °C for short term. Freeze and store at -80 °C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.  
Aliquot the product according to your needs to avoid freeze/thaw cycles.

**ALIAS/SYNONYMS** Teplizumab,MGA-031/MGA031;PRV-031,CD3,anti-CD3

