



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Anti-CLDN18.2 Biosimilar - Anti-Claudin18.2 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1624
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	Mouse IgG1
<b>SOURCE</b>	

### PRODUCT INFORMATION

---

#### ORIGIN SPECIES

**EXPRESSION SYSTEM** Mammalian cells

**BUFFER** PBS buffer pH7.5

**STORAGE CONDITIONS** Store at 4 °C for short term. Freeze and store at -80 °C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.  
Aliquot the product according to your needs to avoid freeze/thaw cycles.

**ALIAS/SYNONYMS** Anti-CLDN18.2,0,CLDN18.2 Biosimilar - Anti-Claudin18.2,anti-CLDN18.2 Biosimilar -  
Anti-Claudin18.2



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**