



Anti NS1 of dengue virus Biosimilar - Anti-Dengue virus mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1627
CLONALITY	Monoclonal Antibody
ISOTYPE	Human IgG1
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer pH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Anti NS1 of dengue virus,0,Dengue virus,anti-Dengue virus

