

Talizumab Biosimilar - Anti-IgE mAb - Research Grade

OVERVIEW

| | |
|-------------------|---------------------|
| PRODUCT REFERENCE | PX-TA1631 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG4-kappa-lambda |
| SOURCE | CAS 380610-22-0 |

PRODUCT INFORMATION

| | |
|--------------------|--|
| ORIGIN SPECIES | Humanized |
| EXPRESSION SYSTEM | Mammalian cells |
| BUFFER | PBS buffer PH7.5 |
| STORAGE CONDITIONS | <p>Store at 4 °C for short term. Freeze and store at -80 °C for long term.</p> <p>Addition of 20 to 40 % glycerol may improve cryoprotection.</p> <p>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.</p> <p>Aliquot the product according to your needs to avoid freeze/thaw cycles.</p> |
| ALIAS/SYNONYMS | Talizumab,C21/AL-90,HU-901,Hu901,TNX-901,IgE,anti-IgE |

