



Acapatamab Biosimilar - Anti-FOLH1;CD3E mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1639
CLONALITY	Monoclonal Antibody
ISOTYPE	Bispecific scFv;Kappa;Lambda
SOURCE	CAS 2314491-93-3

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Acapatamab,ACAPATAMAB,IMMUNOGLOBULIN SCFV-SCFV-SCFC, ANTI-(HOMO SAPIENS FOLH1 (FOLATE HYDROLASE, PROSTATE SPECIFIC MEMBRANE ANTIGEN, PSMA)) AND ANTI-(HOMO SAPIENS CD3E (CD3 EPSILON, LEU-4)), MONOCLONAL ANTIBODY SINGLE CHAIN (SCFV) ² -SCFC, BISPECIFIC,FOLH1;

