



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Barecetamab Biosimilar - Anti-ERBB3 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1646
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1,Lambda
<b>SOURCE</b>	CAS 2275727-74-5

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Homo sapiens
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Barecetamab,IMMUNOGLOBULIN G1, ANTI-(HUMAN EPIDERMAL GROWTH FACTOR RECEPTOR HER3) (HUMAN MONOCLONAL ISU104 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL ISU104 .LAMBDA.-CHAIN, DIMERISU-104,ERBB3,anti-ERBB3



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**