



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Cadonilimab Biosimilar - Anti-PDCD1;CTLA4 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1648
CLONALITY	Monoclonal Antibody
ISOTYPE	Bispecific ,IgG1;Kappa;Lambda
SOURCE	CAS 2394841-59-7

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Cadonilimab,ANTI-PD1/ANTI-CTLA4 BISPECIFIC ANTIBODY AK104, PD-1 X CTLA-4 DUAL CHECKPOINT INHIBITOR AK104,PDCD1;CTLA4,anti-PDCD1;CTLA4



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.