



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Coprelotamab Biosimilar - Anti-ERBB2 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1650
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1,Kappa
<b>SOURCE</b>	CAS 2348469-84-9

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Humanized
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Coprelotamab,IMMUNOGLOBULIN G1, ANTI-(HUMAN TYROSINE KINASE RECEPTOR HER2) (HUMAN-MUS MUSCULUS MONOCLONAL GB221 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN-MUS MUSCULUS MONOCLONAL GB221 .KAPPA.-CHAIN, DIMER, IMMUNOGLOBULIN G1-KAPPA, ANTI-(HOMO SAPIENSERBB2 (EPI



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**