



Dafsolimab Biosimilar - Anti-CD3E mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1652
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2b,Kappa
SOURCE	CAS 2367001-71-4

PRODUCT INFORMATION

ORIGIN SPECIES	Mus musculus
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Dafsolimab,DAFSOLIMAB [INN],IMMUNOGLOBULIN G2B, ANTI-(HUMAN CD3 ANTIGEN .EPSILON..DELTA. COMPLEX) (MUS MUSCULUS MONOCLONAL SPV-T3A .GAMMA.2B-CHAIN), DISULFIDE WITH MUS MUSCULUS MONOCLONAL SPV-T3A .KAPPA.-CHAIN, DIMER, IMMUNOGLOBULIN G2B-KAPPA, ANTI-(CD3E

