



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Daxdilimab Biosimilar - Anti-LILRA4 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1654
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Lambda
SOURCE	CAS 2245966-28-1

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Daxdilimab,DAXDILIMAB, IMMUNOGLOBULIN G1, ANTI-(HUMAN LEUKOCYTE IG-LIKE RECEPTOR A4) (HUMAN MONOCLONAL VIB7734 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL VIB7734 .LAMBDA.-CHAIN, DIMER, MEDI7734, VIB7734, MEDI-7734, VIB-7734,LILRA4,anti-LILRA4



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.