



## Demupitamab Biosimilar - Anti-EGFR mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1655
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1,Kappa
<b>SOURCE</b>	CAS 2249926-74-5

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Homo sapiens
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Demupitamab,IMMUNOGLOBULIN G1-KAPPA, ANTI-(HOMO SAPIENS EGFR (EPIDERMAL GROWTH FACTOR RECEPTOR, RECEPTOR TYROSINE-PROTEIN KINASE ERBB-1, ERBB1, HER1, HER-1, ERBB)), HOMO SAPIENS MONOCLONAL ANTIBODYIMMUNOGLOBULIN IGG1, ANTI-(HUMAN EPIDERMAL GROWTH FACTOR R

