



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Docaravimab Biosimilar - Anti-Rabies Virus Strain ERA GP Ectodomain Epitope G-III mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1658
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2b,Kappa
SOURCE	CAS 2247196-23-0

PRODUCT INFORMATION

ORIGIN SPECIES	Mus musculus
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Docaravimab,DOCARAVIMAB [INN],DOCARAVIMAB [WHO-DD],IMMUNOGLOBULIN G2, ANTI-(RABIES VIRUS STRAIN EVELYN-ROCKITNIKI-ABELSETH GLYCOPROTEIN ECTODOMAIN EPITOPE G-III) (MUS MUSCULUS MONOCLONAL .GAMMA.2-CHAIN), DISULFIDE WITH MUS MUSCULUS MONOCLONAL .KAPPA.- CHAI



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.