



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Ebdarokimab Biosimilar - Anti-IL12B mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1660
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Kappa
SOURCE	CAS 2393651-11-9

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Ebdarokimab,AK-101, EBDAROKIMAB, EBDAROKIMAB, IMMUNOGLOBULIN G1, ANTI-(HUMAN INTERLEUKIN 12 SUBUNIT .BETA.) (HUMAN MONOCLONAL AK101 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL AK101 .KAPPA.- CHAIN, DIMER,IL12B,anti-IL12B



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.