



Erfonrilimab Biosimilar - Anti-CD274;CTLA4 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1665
CLONALITY	Monoclonal Antibody
ISOTYPE	Bispecific Single Domains (VH-VH'-CH),IgG1na;na
SOURCE	CAS 2367013-69-0

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Erfonrilimab,ERFONRILIMAB, ERFONRILIMAB, KN046,CD274;CTLA4,anti- CD274;CTLA4

