



Gomiliximab Biosimilar - Anti-FCER2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1671
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Kappa
SOURCE	CAS 357613-86-6

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Gomiliximab,GOMILIXIMAB, LUMILIXIMAB, IDEC-152, PRIMATIZED ANTI-CD23, ST-152, IMMUNOGLOBULIN G1, ANTI-(HUMAN IMMUNOGLOBULIN E RECEPTOR TYPE II) (HUMAN-MACACA IRUS MONOCLONAL IDEC-152 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN-MACACA IRUS MONOCLONAL IDEC-152 .K

