



Grisnilimab Biosimilar - Anti-CD7 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1672
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a,Lambda
SOURCE	CAS 2367001-70-3

PRODUCT INFORMATION

ORIGIN SPECIES	Mus musculus
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Grisnilimab,GRISNILIMAB,IMMUNOGLOBULIN G2A-LAMBDA, ANTI-(CD7 (CD7 ANTIGEN (P41),GP40, LEU-9, TP41, TP40)), MUS MUSCULUSMONOCLONAL ANTIBOD (GRISNILIMAB), IMMUNOGLOBULIN G2A, ANTI-(HUMAN CD7 ANTIGEN) (MUS MUSCULUS MONOCLONAL WT1 .GAMMA.2A-CHAIN), DISULFIDE

